

Product datasheet for **TP511037**

Prkd2 (NM_178900) Mouse Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Mouse protein kinase D2 (Prkd2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug |
| Species: | Mouse |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >MR211037 protein sequence Red =Cloning site Green =Tags(s) MAAAPSHPAGLPGSPGPGSPPPPGGLDLQSPPLLPQIPAPGSGVSFHIQIGLTREFVLLPAASELAHVK QLACSIVDQKFPECGFYGLYDKILLFKHDPTSANLLQLVRSADIQEGDLVEVVLSASATFEDFQIRPHA LTVHSYRAPAFCDHCGEMLFGLVRQGLKCDGCGLNHYHKCAFSIPNNCSGARKRRLSSTSLASGHSVRLG SSESLPCTAEELSRSTTDLLPRRPPSSSSSSSSSSSYFTGRPIELDKMLMSKVVPHTFLIHSYTRPTVCQ ACKKLLKGLFRQGLQCKDCKFNCHKRCATRPNDCLGEALINGDVPMEEAADYSEADKSSISDELED SGV IPGSHSESALHASEEEEEEGEGHKAQSSLGFIPLMRVVQSVRHTTRKSSTTLREGVWVHYSNKDTLRKRHYW RLDCKCITLFQNNNTNRYEIPLEILAVEPAQNFSLVPPGTNPHCFEITANVTYFVGETPGGAPGGP SGQGTEAVRGWETAIRQALMPVILQDAPSAGHTPHRQASLSISVSNSQIQENVDIATVYQIFPDEVLGS GQFGVYGGKHKRTGRDVAVKVIDKLRFP TKQESQLRNEVAILQSLRHPGIVNLECMFETPEKV FVVM EK LHGDMLEMILSSEKGR LPERLT KFLITQILVALRHLHFKNIVHCDLKPENVLLASAD PFPQVKLCDFGFA RIIGKSFRRSVVGTPAYLAPEVLLNQGYNRSLDMWSVGVMYVSLSGTFPFNEDEDINDQIQNAAFMYP ASPWSHISGAIDLINLLQVKMRKRYSDKSLSHPW LQEYQTWLDLRELEGKMGERYITHESDDARWD Q FVAERHGTPAEGDLGGACLPQDHEMQGLAERISIL TRTRPLEQKLISEEDLAANDILDYKDDDDKV |
| Tag: | C-MYC/DDK |
| Predicted MW: | 97 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |


[View online »](#)

| | |
|----------------------|--|
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C after receiving vials. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_849231 |
| Locus ID: | 101540 |
| UniProt ID: | Q8BZ03 |
| RefSeq Size: | 3625 |
| Cytogenetics: | 7 A2 |
| RefSeq ORF: | 2625 |
| Synonyms: | AI325941; PKD2 |

Summary:

Serine/threonine-protein kinase that converts transient diacylglycerol (DAG) signals into prolonged physiological effects downstream of PKC, and is involved in the regulation of cell proliferation via MAPK1/3 (ERK1/2) signaling, oxidative stress-induced NF-kappa-B activation, inhibition of HDAC7 transcriptional repression, signaling downstream of T-cell antigen receptor (TCR) and cytokine production, and plays a role in Golgi membrane trafficking, angiogenesis, secretory granule release and cell adhesion (PubMed:17226786, PubMed:20819079). May potentiate mitogenesis induced by the neuropeptide bombesin by mediating an increase in the duration of MAPK1/3 (ERK1/2) signaling, which leads to accumulation of immediate-early gene products including FOS that stimulate cell cycle progression (PubMed:17226786). In response to oxidative stress, is phosphorylated at Tyr-438 and Tyr-718 by ABL1, which leads to the activation of PRKD2 without increasing its catalytic activity, and mediates activation of NF-kappa-B (By similarity). In response to the activation of the gastrin receptor CCKBR, is phosphorylated at Ser-244 by CSNK1D and CSNK1E, translocates to the nucleus, phosphorylates HDAC7, leading to nuclear export of HDAC7 and inhibition of HDAC7 transcriptional repression of NR4A1/NUR77 (By similarity). Upon TCR stimulation, is activated independently of ZAP70, translocates from the cytoplasm to the nucleus and is required for interleukin-2 (IL2) promoter up-regulation. During adaptive immune responses, is required in peripheral T-lymphocytes for the production of the effector cytokines IL2 and IFNG after TCR engagement and for optimal induction of antibody responses to antigens (PubMed:20819079). In epithelial cells stimulated with lysophosphatidic acid (LPA), is activated through a PKC-dependent pathway and mediates LPA-stimulated interleukin-8 (IL8) secretion via a NF-kappa-B-dependent pathway (By similarity). During TCR-induced T-cell activation, interacts with and is activated by the tyrosine kinase LCK, which results in the activation of the NFAT transcription factors (By similarity). In the trans-Golgi network (TGN), regulates the fission of transport vesicles that are on their way to the plasma membrane and in polarized cells is involved in the transport of proteins from the TGN to the basolateral membrane (By similarity). Plays an important role in endothelial cell proliferation and migration prior to angiogenesis, partly through modulation of the expression of KDR/VEGFR2 and FGFR1, two key growth factor receptors involved in angiogenesis (By similarity). In secretory pathway, is required for the release of chromogranin-A (CHGA)-containing secretory granules from the TGN (By similarity). Downstream of PRKCA, plays important roles in angiotensin-2-induced monocyte adhesion to endothelial cells (By similarity).[UniProtKB/Swiss-Prot Function]