

Product datasheet for TP511095

Wdr36 (NM_001110015) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse WD repeat domain 36 (Wdr36), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR211095 representing NM_001110015 Red =Cloning site Green =Tags(s)

MAEMESAVEGRTASVLFAGFRALGLFSNEVPHVVRYSALKRRFYVTTTCVGKSFHTYDVQKLSLVAVSNSV
PQDICCMAADGRLVFAAYGNVFSAFARNKEIVHTFKGHKAEIHLLQPGDHSVSDTDSVLIIWHIYSEE
EYLQLTFDKSVFKISTILHPSTYLNKVLGSEQGSLQLWNIKSNKLLYTFPGWKVGVTALQQAPAVDVVA
IGLVSGQVVIHNIKYDETLMKFHQDWGPITSISFRTDGHPIMAAGSPCGHIGLWDLEDKKLINQMRNAHS
TAISGLTFLHREPLLVTNGADNALRIWIFDGPAGEGRLLRFRMGHSASLTIRIYYGQNGQQILSASQDGT
LQSFSTVHEKFNKSLGHGLVNKKIVKRKGLQNTMSVRLPPITQFAAEARQSDWDGIIACHQGKRSCSTW
NYQRSTIGAYFLKPRGVKTNSNLTAVDITSCGNFAIIGLSSGAVDVYNMQSGIHRGNFGDDKAHTGSV
RGVAVDGLNQLVVTAGSERLLKFWNFKSKVLIHSLGLDSSPNMMLLHRDSGILGLAMDDFSIAVLDIETR
KIVREFSGHHGQINDMTFSPDGRWLISAAMDCSVRTWDLPSGCLIDCFLLDSAPLNVMTMSPTGDFLATSH
VDHLGIYLSNISLYSVSLRPLPPDYVPSIVMLPGTCQTQGMEDLEEKTEPSDEMIYESPEQLSEQLV
TLLLPESRWKNLLNLDVIKKKNKPKEPPKVPQSAPFFIPTVPGLVPRFAVPEPSSDPQQSKVNLGILA
QKSNFYLKLEEGLLNNQYEAALNLLKELGPSGIETELRNLSRDDGGGSVEVMRSFLSMIGMMLDRKRDFEL
AQAYLALFLKLHLRTPSEPALLEELVKLSSQVEKDWITHLQSLFNQSMCVLNYIKSAFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	99.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001103485</u>
Locus ID:	225348
UniProt ID:	<u>Q3TAQ9</u>
RefSeq Size:	3281
Cytogenetics:	18 B1
RefSeq ORF:	2697
Synonyms:	5730444A13Rik; Ta-wdrp