

Product datasheet for TP511239

Haus6 (NM_173400) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse HAUS augmin-like complex, subunit 6 (Haus6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR211239 representing NM_173400 Red =Cloning site Green =Tags(s)

MSSARDSAFEKQHLWMYLQALGLDPSSSITFGGKMVPHAHLGENMFDKLNRFDAFHIVSYFLFKTLDEAL
 A
 KEVFRDCWPPFDQKLDMEFRKHCCEWLKEISAECGSSFPQVVGSLLMSPGGPKFIHLMYHFARYVAIKYI
 KTKSNNSLHFAETFNVPQDMHKCLARSHVARNRFLQILQREHYVMQKYQENVNLSVKQVRNARSECM
 SL
 QNQIKRMEPYDEKSNTQEKIQKVRSLWASVNETLMVVEKERDVISSLFGFVNQYALDGAHVAVNIPRLLL
 DRLEEQICRLQIGNVYEAGKLNLLTIHLLNEVLKIMKYEHCAGQAGLTINLHYLEKETKFQRRERLSHM
 KHMHRHKVKETVTAVRQSIADKEGKWHMKWKGLGLSPFRLIKDRPAAVDLLPPMSPLTFDPVSEEVYAKS
 ILFKYPASLPDSHKEPKQESGSRREVDTPGSVCDGANSPVPFPLQCAASSDTNSSTLLEKDPNMRTPIEK
 TLTYKTIADSEEDSPLSDTAKNSQSSVSRGTLPVKRSEPFQKEQARLVDEVARVVLSDSPQLSERKEVE
 LEELIDTLISNPFLTRKQIPRTPENLITDIRSSWRKAVETEDNRSSDGILVDASAREVSSSSSTLHNQR
 EFSLAIFSPDFDRSYFPEEKAVSDQLKGLPQKHLVTSQIREPPTEDVSDLVNKKMIWKQDLECATTRLSE
 TSHMDVFSPPVVGSGIDVMSRSKEDCVPVPPDHSQASCIESTHKTLLWDSFQRLSGIGILHETLPEVGRL
 SLNSSTSSETSFNLEPKSYVHSDAFADDVGKGQAPELDSNLQARLTSYEALKKSLSKIREESYLSNPKT
 LEQHKVELSPVAKAVQADDAHTVWDTQDLFIDYAKPSSSLVPLDGKTQTVSLTKAFLVEQKLKTTVPCSLG
 EFLPSLKEEEISSKSLEARDSLI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	106.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775576</u>
Locus ID:	230376
UniProt ID:	<u>Q6NV99</u>
RefSeq Size:	3719
Cytogenetics:	4 40.69 cM
RefSeq ORF:	2799
Synonyms:	6230416J20Rik; 9830144E06; D4Ert27e; mKIAA1574; RP23-95O1.3