

Product datasheet for **TP511703**

Dhx34 (NM_027883) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse DEAH (Asp-Glu-Ala-His) box polypeptide 34 (Dhx34), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR211703 protein sequence Red =Cloning site Green =Tags(s)
	MPPPRTREGRGHRDRDHHRAPREEEAPEKWDWNCPETRCLLEDVFFRDEDYIRRGSEECQKFWAFFERL Q RFQHLKTSQKKEKDPGMPKHGIAALADLPLTYDPYRINLSILSPDTRGRHGPGRGLPPERVSEFRALL HYLDFQKQAFGRVLKLRERAAPIAQYGNRILQTLKEHQVWVAGDTGCGKSTQVPQYLLAAGFSHVA CTQPRRIACISLAKRVGFESLSQYGSQVGYQIRFESTRSAATKIVFLTVGLLRQIQREPSLPQYQVLIV DEVHERHLHNDLFLGVLQRLLPQRPDLKVLMSATINISLFSSYFSHAPVWQVPGRLPFITVYQVQPEAD QTASKSEKLDPRPFLRVLEAINNKYPPEERGDLLVFLSGMAEITTVLDAAQAYASLTQRWVVLPLHSALS VSDQDKVFDVAPAGVRKCILSTNIAETSVTIDGIRFVVDVSGVKEMSYDPQAKLQRLQEFWISQASAEQR KGRAGRTGPGVCYRLYAESDYDAFAPYPVPEIRRVLDALVLQMKSMSVGDPRTFPIEPPPPASVETAI LYLQEQGALDSSEALTPIGSLLAQLPVDVWIGKMLILGSMFSLAEPVLTIAAALSQSPFTRSAQSNLDC ATARRPLESDQGDPTLNFVNAWVQKSERSGNSRKWCRRRGVEEHRLYEMANLRRQFKELLEDHGLL S GAQVVAPGDSYSRLQQRERRALHQLKRQHEEGGRRRKVLRLQEDGCSSDEEDREGSTSQRADSVDIQ D VKFKLRHNLEQLQAAASSAQDLTRDQLALLKLVLGRGLYPQLAVPDAFNSGRKDSQIFHTQAKQGTVL H PTCVFANSPEVLHTQGQEASGQEGSQDGRDQMSCKHQLLAFVSLLETNKPVLVNCVRIPALQSLLLFSRS IDTNGDCSRLVADGWLELQLADSESAVRLLATSLRLRVHWESALDRQLARQAQRRKLEQEEDMGSPAVSP QEVAALSRELLQFMVAKVPYRLRRLTGLEAQNLVGPQTITTAPSLPGLFGNSTLSPHPTKGGYAVSDYL TYNCLTSDTDLYSDCLRSFWTCPHCGLHMPFTPLERIAHENTCPEAPGDDPGSEEAAPAPPQKTSALQRP YHCQVCGQDFLFTPTVLRHRRQHV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK



[View online »](#)

Predicted MW:	128.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_082159
Locus ID:	71723
UniProt ID:	Q9DBV3
RefSeq Size:	3801
Cytogenetics:	7 A2
RefSeq ORF:	3435
Synonyms:	1200013B07Rik; 1810012L18Rik; Ddx34; mKIAA0134
Summary:	Probable ATP-binding RNA helicase.[UniProtKB/Swiss-Prot Function]