

Product datasheet for TP511771

Heatr6 (NM_145432) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse HEAT repeat containing 6 (Heatr6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR211771 protein sequence Red =Cloning site Green =Tags(s)
	MAAVQVAASLPCGQPREAPRELSPEQDDGFRRLSARLRALQPDDSTVSRMEIHLLFDQLISENYSEGGGV APEDVSALLVRACQLVPLNQNLVSKVCQLIHRLNRLQVVDEPNLDFLLTYTISALQQCSSWTHMEIL QALAALVYCNGSKCQKHLPDLLGKSGLLMKLSDSLHSDPEVRRAAVHCMANLCLSVPGQPYLEEPYQHV C FQAFLTILQSPKSSDMDDITFCMLLQNALKGIQSLNNGGKMRLTQTEHLGALLAVLKKAMFHGLPGLNIE MPAVLYPTPLPQYDGRSPVKPQQPESTTARSSANKKKKHVKTKKTQQGEKAEIIIIIEPCGEVEAAPGLS MDQANVSGESAWCSSPWGSGQLPVDGGRAAGREQASLPFASSSWKRVSSSESDDSDAEGGMQGMRS YQA KVRQGALACFLSTIKSIEKKVLYGYWSAFVPDTPELGSPQSVSLMTLTLKDPSPKTRACALQVLSAILEG SKQFLSVAEDTSDHKRAFTPFSVTIASSIRELHRCLLLALVAESSQTLTQIIKCLANLVSNAPYNRLKL SLLTKVWNHIKPYIRHKDVNVRVSSLTLLGAIVSTHAPLPEVQLLLQQPCSSGLSSSSSATPHLSTPDGW KTFPAGSSLEEASLSSPKGSSEPCWLIRLCISTVVLPKEDSCSGSDAGSALGSTYEPSPMRLEALQVLAH LARGYFSTAQLYLMELGEVICKCMCEANPSIQLHGVKLEELGAGLIQQYKPDSNIAPEQRVPVHLVVVF WTAMLSGPLPRALQSAEHPTLQASACDALSSILPEAFSSLPNDKQILCITVLLGLNDSKNHLVKAATSRA LGVYVLPCLRQDVIFVADTANAILMSLQDKSLNVRKAAWWSLGNLTDTLIVNMDTPDPSFQDEFSGLLL LKMLQCAIQASTDKDKVKSNAVRALGNLLHFLQPSHVERPRFAEIIIEESIQUALISTVNEAAMKVRWNAC YAMGNVFNKPNALPLGTAPWTSQAYKALTSVMSCKNFKVRIRSAALSVPSKRGQYGSLEQFIQIWSALV TALQKSEDTTDFLEFKYCASLRTHICQALLHLLSLASASDLPGIQETLKLHGDMIRSYILQFLKSGAEGD DPGAVHTPQERVQMVRVALKHIHRVQTLAGDTARGAIMGFLEDILTVHCDSSGEQVMLRGPSPDQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	128.9 kDa



[View online »](#)

Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_663407
Locus ID:	217026
UniProt ID:	Q6P1G0
RefSeq Size:	5677
Cytogenetics:	11 C
RefSeq ORF:	3552
Synonyms:	2700008B19Rik; BB135312; BB249351