

Product datasheet for TP512514

Rbbp8nl (NM_173031) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RBBP8 N-terminal like (Rbbp8nl), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR212514 representing NM_173031 Red =Cloning site Green =Tags(s) MDSFMESLNRLKEAHEKEVLGLQNKLELNSERCRDAQRVEELFAKNHQLREQQKALKENLRVLENRLR A GLCDRCMVTQELARKKQLELENAHLQSLQHLCILTNE MNGLREENKILKEEVKRLRSLGDKTAPQAWEGT SEPPSPLPLPSPSSWKGTENPPGAPEEAEEEQSGTDKISSHKTSPVSRVSPAANLPEPRALDMTPQCIS NQLHSTVAVVRPGSRACPPDSGSADGTPPPSTRSTPPSPTYEHGLPVDSFLRASRPSAIAYETLKHSLQT DRLLCLNRQLSLRLQSPHSSPLASAVAPHSPPLRGLKAREAEAWEEP GALVGMQDPRLEGALHLLL VQQ H LRARARASTARLQVPSAPEEMPSSPPAGSDSEGSDSEAPKALSTEAQPDGWHPQTTGQGSSQRKETHVA TQDCPPDKPLDLSDRGRCDISKSTGQPLPLSPKIVYTPSPQLPTLSRPLVPSSHTLSNVSTETRAQESE EYSTPKDTSHPPLPGVCMGLTSPGRTAEEAGGRPQPGPHLQRPDTGGTTEPKKVRAQRLALEQLEESDTS TKAGLISEASAEPSPMPGVGHAEDHQSSQQQKRKRASDLQDKGYLREGKAGPGM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	66.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_766619</u>
Locus ID:	271887
UniProt ID:	<u>A2ABX0</u>
RefSeq Size:	2319
Cytogenetics:	2 H4
RefSeq ORF:	1842
Synonyms:	4732416F18