

Product datasheet for TP513021

Saxo1 (NM_001081096) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse stabilizer of axonemal microtubules 1 (Saxo1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR213021 representing NM_001081096 Red =Cloning site Green =Tags(s) MGKRTCICNLCSGRHRCPLPTKIYDKTEKTRLLSEYTENYPIYQSYLPRNSFKPEWCYRKPSAPMEGL TTCRRDFGLHKMLPVTFHQPDTYVPSQETMDLLTTYKHDFNYLPTCPVGLIKPRDSKFPNGDKIMQYLPT YKVDYLPWNQPRRELLRPPHKYRPESSEKFENRTTHQDDYTMKGLVTTVSCPPSKSKIYNIPLDLTNYK MSYVPHPEKRFVKEAEKYIPCDIPFENLTTHKESYRGLLGEPAKSTKPPAKFPVHDAPFTNTTEVQEKF QAWPTPQMVPKAPVVYPPEEKMDLLTTVQTHYKYLKGSPAITCRPVPSIKKSKCFESSTTTKDDFKQWA AVNTKPIRPTLHLSLPAEPLDCQTTTKTCFVAHPPIITENYKPPWVGLRRNIPVEGQTTYSINFTPKDLG KCPASYMEPPGYIFEEVDMSMGHKIYRSISQTSRQNSSISICDTEKPGQKELEVSA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	55.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001074565</u>
Locus ID:	75811


[View online »](#)

UniProt ID: [B1AXP3](#)

RefSeq Size: 1835

Cytogenetics: 4 C4

RefSeq ORF: 1428

Synonyms: 4930500O09Rik; Fam154a