

Product datasheet for TP513270

Cyp4f39 (NM_177307) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse cytochrome P450, family 4, subfamily f, polypeptide 39 (Cyp4f39), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR213270 protein sequence Red =Cloning site Green =Tags(s)

MLPLTDYLLQLLGLEKTAFRVYVWSALLLVLFFFRLLVRAFKLFSDFRITCRKLSCFPEPPGRHWLLG
HMSMYLPNEKGLQNEKKVLDTMHHIILAWVGPFLPLLVLVHPDYIKPVLGASAAIAPKDEFFYSFLKPWL
GDGLLISKGNKWSRHRLLTPAFHFDILKPYMKIFNQCTNIMHAKWRRHLAEGSVTSFDMFEHISLMTLD
SLQKCVFSYNSDCQERMSDYISSIIELSALVRRQYRLHHYLDYFMYLTDGRRFRQACDVTVHNFTTEVI
QERRRALRQQGAEAWLKAKQGKTLDFIDVLLAKDEEGKELSDDIRAEADTFMFEGHDTTSSGLSWALF
NLAKYPEYQEKREEIQEVMKGRELEELDWDLLTQLPFTTMCIKESLRQFPPVTLSRRCTEDIKLPDGR
VIPKGIICLVSIYGTHHNPVWPDSKVYNPYRFDPTPQQRSPALFVPFSAGPRNCIGQSFAMAEMRVV
ALTLRFRLSVDRTHKVRKPELILRTENGLWLNVEPLPSRA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	62 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_796281</u>



[View online »](#)

Locus ID: 320997

UniProt ID: [Q8BGU0](#)

RefSeq Size: 2618

Cytogenetics: 17 B1

RefSeq ORF: 1599

Synonyms: 4732474A20Rik