

Product datasheet for **TP513445**

Fyb2 (NM_001162980) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse FYN binding protein 2 (Fyb2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR213445 protein sequence Red =Cloning site Green =Tags(s)
	MEARPPERGRSHSPGGAGEEPRLQLLRDRSDCKLAGGGESIVGPRIQLWEAGIVSQSLPRRTTMEEEGI RNFKELRAKFQKFNAAPLPGPMRFPAGVSRKHDRGSTQPAQDLANRCLSSHHHRTPSHSSTGVSQPLK N QTLKSAQGDELQKTSSSSGTPEKSTVCPERDFQKAAVPLDVTKSRAKTSNEEKGMTLSSFRYKLWNWEKV SSQNGEMSPAPLLTNGGIKTFHFEGQKTMGLAQDKPEKSLKATGAQTLPPQSHSMAQRKSPVTSKASSV S LPPHSRKSTKSPATEGSHRSSQCQPVYECELDSLVPKPESRHCQLPKTKPLPSIETLGPPPPKPSKPPF VNLYAFHKLPATVTKPKEETMKEGPLSPDSAELEEAHNYDTTISYLRHSGNSINLCAEGESTEATYEIE IEELQKPWRSFFLPELSPRPKEEENTMEEKESWESEPLEPRKELHPSRPPKVVYKETPGKTQMAGVHED RRSVPAGNQEAMTDIMQNRLFPEVDVTLTRHSQDKSGYVEALEVTKETPSPSTIRSSSSSEKTYDAVECSR EDVRKWDFSSSFTSDSEENCEEMYEDIYKAKNDDPKTEVAGRTALRKLQQIFRKENVVFRMKKTKSKEIV SNGFSVSLPDLGPRSQDDSDGIIYDDVDTREKESNDEDKVKTWTKFLIPKGKKWGKGSQGSKSFSPRH FFRTKKQKLEKNRMEKEEKLFRERFQYGKEILVINRAVACASNSRNGMFDLPPIPGEQLEVIDTTEQNLV ICRNSKGKGYVLVEHLDFKHQG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	89.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001156452
Locus ID:	242594
UniProt ID:	E9PZ08
RefSeq Size:	3136
Cytogenetics:	4 C6
RefSeq ORF:	2379
Synonyms:	1700024P16Rik; ARAP; C1orf168; Gm620