

Product datasheet for TP513489

9930104L06Rik (NM_177573) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RIKEN cDNA 9930104L06 gene (9930104L06Rik), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR213489 protein sequence Red =Cloning site Green =Tags(s) MAHDQPLLWQEALRKCFPWEEQQNLWQSTLQDCSPLLSSLSNLAEQLQAAQSLRFEDVPALRPFPDL Q ERLRRKQLEAGDVWLDKLAERLATLLKVRNTINSHVEQVFQAYEQHAAVLDTVLRPSVWSPSVADMLE WLQDIDRHYGSSYLKRKYLSSIHWDLASIQALPKAWDQISENECQTLVSDVLVSVSFFLEEPGGCAAS GDLEHHS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	24.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_808241
Locus ID:	194268
UniProt ID:	Q499E6


[View online »](#)

RefSeq Size:	2942
Cytogenetics:	4 D2.2
RefSeq ORF:	651
Synonyms:	C1orf109; MGC107426