

Product datasheet for **TP513510**

Fam163a (NM_177838) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Mouse family with sequence similarity 163, member A (Fam163a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

Species: Mouse

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >MR213510 protein sequence
Red=Cloning site **Green**=Tags(s)

MTAGTWWITGGILATVILLCCIIVLCYCRLQYYCCKKGTGDGEDAESEEEEEEEHGLSIHPRVPACNACSSH
VLDGRGGLAPLTSESCSQPCGVASHCTTCSFYRTPFYIRTADMVPNGGGGERLSFAPTHYKEGGTPSLKL
AAPQNYPTWPSSGHEAFTNPRAISTDV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-MYC/DDK

Predicted MW: 17.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_808506](#)

Locus ID: 329274

UniProt ID: [Q8CAA5](#), [A9ZNB6](#)

RefSeq Size: 3892

Cytogenetics: 1 G3



[View online »](#)

RefSeq ORF: 507

Synonyms: A230106N23Rik