

Product datasheet for TP513783

Smco3 (NM_001039558) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse single-pass membrane protein with coiled-coil domains 3 (Smco3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR213783 representing NM_001039558 Red=Cloning site Green=Tags(s) MAQSDFLYPQNPRRRQEVNRLHQQLLDCLSDSFQVTNKLTVLNLHLCRLAFIEMKSDGTIKENCDIII QAMTKIQKELQKIDEALKDQLEPTLYRKLQDIKERETEKIAIVQKVISVILGEATSAASAVAVKLVGSSV TTGIISKLVSVLAHIGTSLLGSIGVAVLSLGIDMIIQAILGAVERTQLQAAIKSYEKHLEEFKAASAKYH HAITEVATAVKRQLR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	25.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001034647</u>
Locus ID:	654818
UniProt ID:	<u>Q8BQM7</u>
RefSeq Size:	2052


[View online »](#)

Cytogenetics: 6 G1

RefSeq ORF: 675

Synonyms: C030030A07Rik