

## Product datasheet for TP514127

### Ccdc24 (NM\_001034876) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse coiled-coil domain containing 24 (Ccdc24), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR214127 representing NM_001034876 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MLKALLQEIQSSQVNSDPPSLLAPPPLLRDLIRQELRQLLQGLRLKAISEGRDQTQMWAQYSPRVLRFAL LEDTRCDSTQQELPMRADEPRAGMNNTERVMVTSPVSTCPRLDQVQVVRHLRSLLEEECHMLQNEISDLQHC LEMEQQACQPSKDTQMPTLGEIKDQKAAMERELQVSLGPSCTAAKHRQKPLRSSIP SLRPFPCPLHGYMPTPPSERCPHPQGQTSTRRWGRQLRYSYRKETTSTSVSSAASQAPT  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	30.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001030048</a></u>
Locus ID:	381546
UniProt ID:	<u><a href="#">Q3V0I8</a></u>


[View online »](#)

RefSeq Size:	1295
Cytogenetics:	4 D2.1
RefSeq ORF:	804
Synonyms:	Gm1028