

Product datasheet for **TP514141**

Scoc (NM_001039137) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse short coiled-coil protein (Scoc), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from Mouse cDNA ORF Clone, MR214141, encoding Mouse full-length Scoc.
Tag:	C-MYC/DDK
Predicted MW:	14.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001034226
Locus ID:	56367
UniProt ID:	Q78YZ6
RefSeq Size:	1835
Cytogenetics:	8 C2
RefSeq ORF:	375
Synonyms:	0610011C06Rik; 2810002J11Rik; AA415030; SCOCO
Summary:	Positive regulator of amino acid starvation-induced autophagy.[UniProtKB/Swiss-Prot Function]



[View online »](#)