

Product datasheet for TP514181

Upf3a (NM_025924) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse UPF3 regulator of nonsense transcripts homolog A (yeast) (Upf3a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR214181 protein sequence Red =Cloning site Green =Tags(s) MRSEEGAGGLGAALAARGPSWREKLSSSETHCRRESPRKESAAPPAPTLSSEGAGKPREEKRTALSKVVL RRLPPGLTKEQLEEQLRPLPAHDYFEVVAADLSLYPHVYSRAYINFRNPDDILLFRDRFDGYIFIGNKGL EYPAVVEFAPFQKIAKKKKLKKDAKTGSIEDDPEYKQFLESYSLEEEKTSASPETLLGEIAKTRELLAR RTTPLLEYIKNRKLEKQRLREEKREERRRRLEKKRLREEEKRKRREEDRCKRKEAERHKRAEKDVRIKL LKKLETGEEVATEKPKERGEAVGAGDEKQEDRLVASQPETAGERLQDTSKEPRDKRSREKDPETSRCHP DTCRKHSSHSEFARFSRRNEDEAKWGKGFAQDRGKKGQDGSVPMDDVAERSGKEHKEYDSLGPDKDRLG HK AQ TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	48.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_080200</u>


[View online »](#)

Locus ID: 67031

UniProt ID: [Q3ULJ3](#)

RefSeq Size: 1402

Cytogenetics: 8 A1.1

RefSeq ORF: 1266

Synonyms: 2600001C03Rik; 4930546M19Rik; RENT3A; UPF3