

Product datasheet for TP514184

Clhc1 (NM_028504) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse clathrin heavy chain linker domain containing 1 (Clhc1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR214184 representing NM_028504 Red =Cloning site Green =Tags(s) MSVQEVNTHAVLPPIVSRNDKEFLESIQRYITTETKRVGCKEEGPADEYYTIYRNVFDKFSATDRKIIFQ VIDYVSAYKSILTAIKKEYDAFIETIKKGRRTAFYLGKLVLAKEPTALVYHQRRAIQLEAKMRIENN STAIQLQIDQMKQLRMEYDKKEVKLCAPSRQLWKPIPGMTLQDSVNLEALNKHKQYLEDKYIKLKQDMS T MYVPAQKKAELDEEMVLLNRRDIAENLKKDRQFRHQRLQVISHTLTPWMKQNMRI SFQDVMERIRKTK A IYGDNIVDEIFEDDPNKKKEAIVMLHYIERFNDLISLGEYERAA CF AANSPKRILQNTSTMNKFKAIGK IRGKPLPLLLFFEAIFNTSQA FKR PINADLTMEGIKCG LSEERLDLVTHWVTQEKLTFSEKAGDIIFAYG EQHTYHKPRCLALAQIIYNECGLHRKALLCLCKQGQIHEAMEHIQQSKDINTDDLIQLITACPQIDLIRC LTQERNEKPPFLSFG L AVLHMFSDMKKVGMRLLQEVSKEKDVIEHLVMSDLFCSLEKWQEIANICLQN GFKILFNDIMSILRSQAGVSEIS EDDTTNIMEHVFVW TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	69.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_082780
Locus ID:	73324
UniProt ID:	Q5M6W3
RefSeq Size:	2185
Cytogenetics:	11 A3.3
RefSeq ORF:	1788
Synonyms:	1700034F02Rik