

Product datasheet for **TP514346**

Ankrd34a (NM_001024851) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Mouse ankyrin repeat domain 34A (Ankrd34a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

Species: Mouse

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >MR214346 protein sequence
Red=Cloning site **Green**=Tags(s)

MLHTEGHALLRAVGQGKRLRLRLLEGGAYVNEGDAQGETALMAACRARYDDPQNKARMVRYLLEQGA
DP
NIADRLGRTALMHACAGGGGAASLLAHGADPSVRDHAGASALVHALDRGDRETLATLLDACKAKGTE
VIIIITDTSPSGTKKTRQYLNSPPSPGVDPAPAPPSPGVCTSPSEVQLQTAGGGGRGLLSPRAQEEEEKR
DIFEFLPKPPDDPSPSEPLPKPPRHPPKPLKRLNSEPWGLVAPPQVPPAEGRPGLERLAAEFNGLTLT
GRPRLSRRHSTEGPEDPPPWAKEVTGGGPLSRRNTAPEAQESGLPSGLRQKLSRMEPVELDTPGHFCPDS
PESSRLSLERRYSASPLTLPAGSVSSPRQSQESLPGAVSPLSGRRRSPGLLERRGSGTLLLDHISQTR
PGFLPPLNVSPHPIPIRQPGGRAPSLPAPPPSGAPGSPRTKRKLVRHSMQTEQIRLLGGFQSLGGP
GEPGR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-MYC/DDK

Predicted MW: 52.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq:	<u>NP_001020022</u>
Locus ID:	545554
UniProt ID:	<u>B2RW11</u>
RefSeq Size:	2586
Cytogenetics:	3 F2.1
RefSeq ORF:	1485
Synonyms:	Ankrd34; AW060394