

Product datasheet for TP514459

Zmynd12 (NM_001014900) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse zinc finger, MYND domain containing 12 (Zmynd12), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR214459 representing NM_001014900 Red =Cloning site Green =Tags(s) MNVIYPLAVPKGRRLSCEVCEAPAERMCTACTVTYYCGPVHQKADWDSIHAKICQLLIPLRTTMPFYNSE EERQHGLQQLQRRQKHLIEFCYTVAQKYLFEGRHEDAVPAALHSLRFRMNVHGLSSVELVPAYLLAEAS LGLGRIVQAEEYLSQAQWTVLKSTDCSNATHSLLHRNLGLLYIAKKNLDEARYHLANDIYFASCAYGTED IRTSGGYFHLANIFYGLKKMELADTLYTKVSEIWYTYLNSHYQALLEARSQQTDLLGKQFVNDTGLDEAQ EAEAIRILTSTLIQESTSTKAPQKSILILKTL SILYYLMLECSKAKEYAMKALSLAETVSEEEQRFIQ ELLSVISAEEKQP TR TRPLEQKLI SEEDLA ANDILDYK DDDDKV
Tag:	C-MYC/DDK
Predicted MW:	41.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001014900
Locus ID:	332934
UniProt ID:	A2BGJ5


[View online »](#)

RefSeq Size:	1505
Cytogenetics:	4 D2.1
RefSeq ORF:	1089
Synonyms:	AI414736