

Product datasheet for **TP514918**

Ccdc178 (NM_027616) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse coiled coil domain containing 178 (Ccdc178), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR214918 representing NM_027616 Red =Cloning site Green =Tags(s)

MPENEKEPAQPTTNE DALDTGKTCLQIKELKDKERALMLAENARMQGLSEFFFSRHEDKRLYESKSTNTE
DVNKAIFYSPSRHSCTLVNIPKPCVNKMISHIEDVESKIQEHLTQFEASFEEWTSVTKDKEAGLDVSA
PEKQVHPEKGKDEKCPCLKRMETLLSEAIHLIKSLETDRAEAEQALKQHKSRRKKISMKIDSWSIWKLO
ELPLAVQKEHENFSKINAELRSYLEDIALKVEQLQERKEKLEKANAKLQVDIDYMASHSVLLEKKRKQEL
GCLKERYHKKFEVMEKFRAIHEELKESVDKCEGAKARLKNMKVENEREIQEELINATSYEKELDKLSVLD
AHYTTSIETVNLIDIEGDEEAMNEVLRETQSTTNELENLKKTVDDLKRLFDQYCWRQRKYENEYLEAFSNF
YSLKKTWDIELSNVSKDAKDLTIVYEAQSEENKRIQSEIQSITDDIEESIKKTAEMEEVHTFVEMKMKN
NNYLKQLYKQAYQVGAVYHLSRHKTEELEDKLADLKRIFKGREELLKKLTRGDIATGIEIQRLYAIET
QFIEMQEFIRRVLYNMALLEVEEQLKELEAEAVRIRYLHRQHSKMLHNIRKRKERVKNVNDATKKKLLK
KSKKSRMELTRTEGKRSHHEEIEIARGQTVLHEKCIELSKEIRIMNLERTNYEERLKKLQEEFFKLQF
DREHVGVDHLMREKQYCEERIFEEERRFRRLIDMRKNTLKNIRKCQDDLLEENLRLAKEYQSAQLIFL
KEKESYFNGYDRLLSLNFSLSDKKKLCLQKRLDQKWQEYFRLMILFNKTKLAKFQGDSQNSIQKILAVQ
EESSLMQHILDFFKSPNSSCGEDD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	103 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_081892</u>
Locus ID:	70950
UniProt ID:	<u>Q8CDV0</u>
RefSeq Size:	3335
Cytogenetics:	18 A2
RefSeq ORF:	2598
Synonyms:	4921528I01Rik