

Product datasheet for TP515483

Atpaf1 (NM_181040) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ATP synthase mitochondrial F1 complex assembly factor 1 (Atpaf1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR215483 representing NM_181040 Red =Cloning site Green =Tags(s) MEQARGPGAADAPEEEEGESQAAMAAVSAAGGACPAVLQVAGLYRGLCAVRSRALGLGFVSPAQLRVF PVRRGSGLPPEGADGSGVSELEANPFYDRYDKIQQLRRSDPAAFESRLEKRSEFRKQPVGHSKQSDFIK CMEQKTDALGKQPVSKGFTKDKTLSSVFNEMVKDKTAAEIKQIWQQYFSAKDTVYAVIPKEKFDLIWNR AQSCPTFLCALPRRDGYEFFVGQWTGTTELHFTALINIQTRGDAAASQLILYHYPELKEEKGIVLMTAEMD STFLNVVEAQCIANQVQLFYATDRKEIYGLVETFNFRPNEFKYMSVIAELEQSGGLGAELKRAQNQDKT TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	38.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_851383
Locus ID:	230649
UniProt ID:	H3BLL2


[View online »](#)

RefSeq Size:	2330
Cytogenetics:	4 D1
RefSeq ORF:	1044
Synonyms:	6330547J17Rik; AI593252; ATP11; ATP11p