

Product datasheet for TP515741

Lingo2 (NM_175516) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse leucine rich repeat and Ig domain containing 2 (Lingo2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR215741 representing NM_175516 Red =Cloning site Green =Tags(s)

MLHTAIPCWQPFLGLAVLLLMGSTIGCPARCECSAQNKSVSCHRRRLAIPGPIETKILDLSKNRLK
SINPEEFISYPLLEEIDLSDNIIANVEPGAFFNLFNLRSLRLKGNRLKLVPLGVFTGLSNLTKLDISENK
IVILLDYMFQDLHNLKSLEVGDNLDVYISHRAFSGLLSLEQLTLEKCNLTAVPTEALSHLRSLIALHLKH
LNINNMVPVYAFKRLFHLKLNLEIDYWPLLDLMPANSYGLNLTSLTNTNLSTVPFLAFKHLVYLTHLNL
SYNPSTIEAGMFSDLIRLQELHIVGAQLRTIEPHSFQGLRFLRVLNVSQNLLLEENVFSSPRALEVL
SINNNPLACDCRLLWLLQRQPNLQFGGQQPMCAGPDTIRERSFKDFHSTALSFYFTCKKPKIREKKLQHL
LVDEGQTVQLECNADGDPQPVISWVTPRRRFITTKSNGRATVLGDGTLEIRFAQDQDSGMVYCIASNAAG
NDTFTASLTVKGFTSDRFLYANRTPMYMTDSNDTVSNGTNANTFSLDLKTLVSTAMGCFTFLGWVLFCE
LLLVWSRGKGGKHKNSIDLEYVPRKNNGAVVEGEVAGPRRFNMKMI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	68.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_780725](#)
Locus ID: 242384
UniProt ID: [Q3URE9](#), [B2RQD4](#)
RefSeq Size: 3669
Cytogenetics: 4 A5
RefSeq ORF: 1818
Synonyms: B230217C06Rik; Lern3; Lrrn6c