

## Product datasheet for **TP515933**

### Rab4a (NM\_009003) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RAB4A, member RAS oncogene family (Rab4a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR215933 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSETYDFLFKFLVIGNAGTGKSCLLHQFIEKKFKDDSNHTIGVEFGSKIINVGGKYVKLQIWDTAGQERF RSVTRSYRGAAGALLVYDITSRETYNALTNWLTARDMLASQNIIVLILCGNKKDLADREVTFLASRFA QENELMFLETSALTGENVEEAFMQCARKILNKIESGELDPERMGSGIQYGDAAALRQLRSPRRTQAPSAQE CGC
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	23.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_033029</a>
Locus ID:	19341
UniProt ID:	<a href="#">P56371</a> , <a href="#">Q3ULK1</a>
RefSeq Size:	1427



[View online »](#)

<b>Cytogenetics:</b>	8 72.19 cM
<b>RefSeq ORF:</b>	642
<b>Synonyms:</b>	AI848268; Rab4
<b>Summary:</b>	Protein transport. Plays a role in vesicular traffic. Mediates VEGFR2 endosomal trafficking to enhance VEGFR2 signaling.[UniProtKB/Swiss-Prot Function]