

Product datasheet for **TP516019**

Megf6 (NM_001162977) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse multiple EGF-like-domains 6 (Megf6), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR216019 representing NM_001162977 Red =Cloning site Green =Tags(s)
	MPVGVEARASWRVVALTLLLLPAVPAASQPLPPRPLQPSMPHVCAEQKLTLVGHRQPCVQAFSRVVPVW R SGCGQQAQWCVGQERRTVYYMSYRQVYATEARTVFRCCPGWSQKPGQEGCLSDVDECANANGGCEGPCC NT VGGFYCRCPPGYQLQGDGKTCQDVDECRSHNGGCQHRCVNTPGSYLCECKPGFRLHTDGRTCLAISSCT L GNGGCQHQCQVQLTVTQHRCQCRPQYQLQEDGRRCVRRSPCADGNGGCMHTCQELRGLAHCGCHPGY QLAA DRKACEDVDECALGLAQCAHGCLNTQGSFKCVCHAGYELGADGRQCYRIEMEIVNSCEAGNGGCSHGCS H TSTGPLCTCPRGYELDEDQKTCIDIDDCANSPCCQQVCANTPGGYECSCFAGYRLNTDGCAGEDVDECAS GHSGCEHHCSNLAGSFQCFCEAGYRLDEDRRGCTPLEESVVDLDGQLPFVRPLPHIAVLGDELPQLFQDD YGAEAAVAELRGEHTLTEKFVCLDHSFGHDCSLTCDDCRNGGTCFPGLDGDCDCPEGWTGIICNETCPPD TFGKNCSSPICQNGGTCDPVSGACRCPPGVSGAHCEDGCPKGFYKHKCRKKCHCANRGRCHRLYGACL C DPGLYGRFCHLACPPWAFGPGCEDCLCEQSHTRSCNSKDGSCSCKAGFQGERCEAECEPGSFGPGCRN R CTCRPGVACDPVSGECRTQCPPGYQGEDCGQECPVGTFGVNCSSGSCSCVGAACHRVTECLCPPGKTGE D CGADCPEGRWGLGCQEICPACEHGASCDPETGTCLCLPGFVGSRCQDACPAGWFGTGCQMRCACANDG HC HPATGRCSAPGWTGLSCQRACDSGHWGPDCIHPCNCSAGHGNCDAVSGLCCLCEAGYEGPQCEQWCR QGY FGPGCEQKCRCEHGATCDHVSACTCPAGWRGSGFCEHACPAFFGLDCGSACNCSAGAPCDAVTGSCIC P



[View online »](#)

AGRWGPHCAQTCPPPLTFGLNCSQICTCFNGASCDPVLGQCHCAPGWMGPTCLQACPAGLYGKNCQHS
 CLC
 RNGGNCDPILGQCTCPEGWTGLACENECLPGHHGAGCRLNCSCLNGGTCDRLTGHCRCPTGWTGDKC
 QSP
 CVSGMFGVHCEEHCACRKGATCHHVTGACLCPPGWRGSHCEQACPRGWFGCAQACRCHCPPGASCHH
 VSG
 ECHCPPGFTGPGCEQACQPGTFGKDCEHPCQCPGETWACHPASGACVCAAGYHGTDCQQRCPSGRYGP
 GC
 EQICKCLNGGTCDPATGACYCPAGFLGADCSLACPQGRFGPSCAHVCTCGQGAACDPVSGTCICPPGKTG
 GHCEGCPQDRFGKGCEHKACRNGGLCHATNGSCSCPLGWMGPHCEHACPAGRYGAACLLECSCQN
 NGS
 CEPTSGACLCGPGFYGQACEDTCPAGFHGSGCQRVCECQQGAPCDPVSGRCLCPAGFHGQFCERGCKPG
 F
 FGDGCLQQCNCPTGVPCDPISGLCLCPPGRAGTTCDLDCRRGRFGPGCALRCDGCGGADCDPISGQCHC
 V
 DSYTGPTCREVPTQLSSIRPAPQHSSSKAMKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	165.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001156449
Locus ID:	230971
UniProt ID:	A0A0A0MQ83
RefSeq Size:	6876
Cytogenetics:	4 E2
RefSeq ORF:	4716
Synonyms:	2600001P17Rik; Egfl3