

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_598329
Locus ID:	81910
UniProt ID:	Q99PL5 , Q8CAB9
RefSeq Size:	2706
Cytogenetics:	2 G1
RefSeq ORF:	2103
Synonyms:	1700087N07Rik; 5730465C04Rik; ES/130; mKIAA1398; mRRp0; mRRp1.8; mRRp10; mRRp15a; mRRp15b
Summary:	Acts as a ribosome receptor and mediates interaction between the ribosome and the endoplasmic reticulum membrane.[UniProtKB/Swiss-Prot Function]