

Product datasheet for **TP516076**

Plcxd1 (NM_207279) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse phosphatidylinositol-specific phospholipase C, X domain containing 1 (Plcxd1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR216076 representing NM_207279 Red=Cloning site Green=Tags(s) MSQSELVPTPLMQGLEPSRVGRVGEVCRSVHKNVLTVP EHTVGTEAAAPMGDQADWMSQLCPQLWDVPL HHLSIPGSHDTMTYCLNRKSRISRASSWLLHLLGRVVPFITGPVVMKWSVTQTL DVTQQLDAGVRYLDLR IAHAPGSTRNLCFVHMMYTKALVEDTLTEIAEWLQSHPREVVILACRNFEGMTCELHDYLAGCIVNIFG DMLCPSGEVPTLRQLWAREQQVIVSYEDEATVSRYDQLWPAIPYWWGNAVKTDVLLR FLETMKGGQGRPD G LFVAGINITENLCYILLHPVDSLEEMTRRSLPLMTEWVCAQQPGQSPQCTNIIAGDFVDADGFVSKVISL NCKLLSP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	40 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_997162</u>


[View online »](#)

Locus ID: 403178

UniProt ID: [G3X9J3](#)

RefSeq Size: 3484

Cytogenetics: 5 F

RefSeq ORF: 1071

Synonyms: A330045H12Rik; Gm576