

Product datasheet for TP516256

Anxa3 (NM_013470) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse annexin A3 (Anxa3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR216256 representing NM_013470 Red =Cloning site Green =Tags(s)
	<p>MASIWVGPRGTIKDYPGFSPVDAEAIKKAIRGLGTDEKTLINILTERSNAQRQLIVKQYQAAAYEQELKD DLKGDLSGFHFEHVMVALVTAPALFDAQQLKKSMTGTGDEDALIEILTRSSRQMKEISQAYYTVYKKS GDDISSETSGDFRKALLTLADGRRDESLKVEHLAKKDAQILYNAGENKWGTDEDKFTVLCRLRFPQLK LTFDEYRNISQKDIEDSIKELSGHFEDLLAIVHCARNTPAFLAERLHQALKGAGTDEFTLNRIMVSRS EIDLLDIRHEFKKHGYSLSAIQSDTSGDYRTVLLKICGEDD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	36.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_038498
Locus ID:	11745
UniProt ID:	O35639



[View online »](#)

RefSeq Size:	1619
Cytogenetics:	5 47.29 cM
RefSeq ORF:	969
Synonyms:	Anx3
Summary:	Inhibitor of phospholipase A2, also possesses anti-coagulant properties.[UniProtKB/Swiss-Prot Function]