

Product datasheet for TP516383

Dhrsx (NM_001033326) Mouse Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Mouse dehydrogenase/reductase (SDR family) X chromosome (Dhrsx), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug |
| Species: | Mouse |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >MR216383 representing NM_001033326 Red=Cloning site Green=Tags(s) |

MSTLRALRAVLVCVYAVGIAVALAQLLRRLRGDFRPPVLPQPGRVAIVTGATAGIGRSTARQLARLGMCV
VVAGNDEHCGQEVVSSIRAEMGSDRAHFLPLDLASLASVRGFARDFRALGLPLHLLVNNAGVMLEPRAET
EDGFERHLGVNFLGHFLLTLLLLPALRASGAEGRSRVTVGSATHYVGTVDMADLHGRHAYSPYAAYAQ
SKLALALFALQLQRILDARGDPVTSNMADPGVVDTELYRHAGWVLRVTKRFLGWL VFKSPEEGAWTLVYA
AAAPELEGVGGRYLRDEAEAEPLGTARDQELQRRLWAEGLRLTGAGGGDSGDVAV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-MYC/DDK |
| Predicted MW: | 36.4 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C after receiving vials. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | NP_001028498 |
| Locus ID: | 236082 |
| RefSeq Size: | 2211 |



[View online »](#)

Cytogenetics: 4
RefSeq ORF: 1005
Synonyms: C80877; D4Ert296e; Dhxsxy; Pscad