

## Product datasheet for **TP516672**

### Rab37 (NM\_021411) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse RAB37, member RAS oncogene family (Rab37), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR216672 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MTGTPGAATAGDGEAPERSPPFSPNYDLTGKVMLLGDSGVGKTCFLIQFKDGAFLSGTFIATVGDIFRNK VVTVDGARVKLQIWDTAGQERFRSVTHAYYRDAQALLLYDITNQSSFDNIRAWLTEIHEYAQRDVMIML LGNKADVSSERVISEDGETLAREYGVPFMETS AKTGMNVELAFLAI AKELKYRAGRQPDEPSFQIRDYV ESQKKRSSCCSFV  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	24.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_067386</a>
Locus ID:	58222
UniProt ID:	<a href="#">Q9JKM7</a> , <a href="#">Q544E8</a>
RefSeq Size:	2193



[View online »](#)

Cytogenetics: 11 E2

RefSeq ORF: 672

Synonyms: B230331O03Rik; B230354I04Rik