

Product datasheet for **TP516795**

Dmrtb1 (NM_019872) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse DMRT-like family B with proline-rich C-terminal, 1 (Dmrtb1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR216795 representing NM_019872 Red =Cloning site Green =Tags(s) MLRAPKCSRCRNHGVLVPVKGHTGKCRWKQCICDKCYLITERQKIMAAQKVLRTQAAEEQVATVGTQGPQ LPPRAPAAAATALSSICPLPRAVPGGVGPGPTATCFLERPPQAPSPGPSTFQLGPSGRPGPSTFQPGPG APGGLRDRSSAWLPQLMPQAPRPELCYPDQHLVPVRVPVPGPVRPVRPLPFADYGHPLRFKSDHVVGAGN PEREPFKQCPACVPVSPYQSFPLSEGQDSSSALGVPPQQRGFRHVSCSPYHRSGLVSEPARDLQPTYCSPP PPPPPPPPPLAPPPQPQQPHFLPPGYLSALHFLPPPPPPSPPSFSLTYDTDKENTNDQDAEAPTEPS QDSPQEASN TR TRPLE QKLISEEDLA ANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	38.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_063925
Locus ID:	56296
UniProt ID:	A2A9I7
RefSeq Size:	1797
Cytogenetics:	4 C7
RefSeq ORF:	1077
Synonyms:	Dmrt6; Prpl3