

Product datasheet for TP517293

Pih1d2 (NM_028300) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse PIH1 domain containing 2 (Pih1d2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR217293 representing NM_028300 Red =Cloning site Green =Tags(s)
	 MTVSSKGLLTHISQFWNMLDDLAENDPERYRNFIQQELKDQKQLCVNPEPQLCIQTKILKPNEKVLFINL CQWERIPAPQSATRPVPSVGRPEDSAEASDAYTIIDVAYNPGVLQAAEKDQGIKDLIRMAMLCIEERL QFTLAHSYHLTSFRLKGSIQRMKESLMGIKDFPDLKAMTRTENTLARIRSSTVSEPNHLPEVLLTKKQA SAKGRCLIEEISSSEIQVEVKKPAYELKVVKDRNEKPLKIELKVELPGIKSVSLCELSVSEVDILIEVSE MYRLCLNLPESINTEMTTAKFVKNKKSALIITMPLA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	35.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_082576
Locus ID:	72614
UniProt ID:	Q8CHR9



[View online »](#)

RefSeq Size: 1268
Cytogenetics: 9 A5.3
RefSeq ORF: 945
Synonyms: 2700059L22Rik