

## Product datasheet for TP517360

### Polr1c (NM\_009085) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse polymerase (RNA) I polypeptide C (Polr1c), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR217360 representing NM_009085 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAAAQAVEEMRTRVWLGEFGVRNVHTTDFPGNYAGYDDAWDQNRFEKNFRVDVWQMDEDTLEFDMVGIDA AIANAFRRILLAEPVPTMAVEKVLVYNNTSIVQDEILAHRLGLIPILADPRLFEYRNQGGEEGTEIDTLQF RLQVRCTRNPNAAKDSSDPNELYVNHKVYTRHMTWVPLGNQADVPEGTIRPVHDDILIAQLRPGQEIDL MMHCVKGIGKDHAKFSPVATASYRLLPAITLLEPVEGEEAAEELSQCFSPGVIEVEEVQGKKVARVANARL DTFSREIFRHEKLLKAVRLARVRDHYIFSVESTGVLPPDVLVSEAIKILMGKCRFLDELDAVEMD</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-MYC/DDK
Predicted MW:	39.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_033111</a>
Locus ID:	20016
UniProt ID:	<a href="#">P52432</a>



[View online »](#)

**RefSeq Size:** 1305

**Cytogenetics:** 17 C

**RefSeq ORF:** 1038

**Synonyms:** 40kDa; AA409007; AA959927; AL024089; Polr1e; RPA40; Rpo1-1

**Summary:** DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Common component of RNA polymerases I and III which synthesize ribosomal RNA precursors and small RNAs, such as 5S rRNA and tRNAs, respectively. RPAC1 is part of the Pol core element with the central large cleft and probably a clamp element that moves to open and close the cleft (By similarity).[UniProtKB/Swiss-Prot Function]