

Product datasheet for TP517381

Txlnb (NM_138628) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse taxilin beta (Txlnb), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR217381 protein sequence Red =Cloning site Green =Tags(s) MEINHPDQLSVEHPTPPGDSSSLNQNGPGKQDGERCSTSGQAPEQEGSLHPEKGAHDVAEELSRQLEDII STYGSAASPRGKESTSETKEQPPNTEAPDNEDVDYEETTEEIDREPTAPEEPAAAKEPVSNEQKLEKKI LKGLGKEANLLMQNLNKLQAPEEKLDLFLFKYTELLDEHRTEQKKLKKLLKQQAQTQREKDLQSEHNR A VLARSKLESCLRELQRHNKTLKEETLQRAREEEERKEITSHFQTTLDIQTQIEQQSERNMKLCQENTE LAEKLKSIIDQYELREEHLDKIFKHRELQKLVDAKLEEAQELMQEAERHRREKEYLLNQAAEWKLQAK VLKEQETVLQAQLTLYSGRFEFQSTLTKSNEVFATFKQEMDKTTKMKKLEKDTATWKARFENCNKALL DMIEEKALRAKEYECFVMKIQRLENLCRALQEERKELYKKIREAKMSEKEDQVQRTSEEEPEPSVSENEE VDAEEANSFQKAVERNATAFTILHHPEFTPDQPTERQLAVNGPQSGSDVTHQHPETARLNHPSLPADSGS PRPPVGAQAVAEGVCEATPAPTASCTPAEALQSQGLPAENTPGPKPHKPEANTSGQAPLSPAQGSLSVV EAKYDISPSESEGDSAVVPGCESREQPPPEVTDIPVGPSTGLPREPDACLNQVD TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	77.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_619534</u>
Locus ID:	378431
UniProt ID:	<u>Q8VBT1</u>
RefSeq Size:	4457
Cytogenetics:	10 A2
RefSeq ORF:	2055
Synonyms:	2310001N14Rik; Mdp77
Summary:	Promotes motor nerve regeneration. May be involved in intracellular vesicle traffic (By similarity).[UniProtKB/Swiss-Prot Function]