

Product datasheet for TP517396

Cenpl (NM_001159930) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse centromere protein L (Cenpl), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR217396 representing NM_001159930 Red =Cloning site Green =Tags(s) MDSCDFRGLTSRRRTSNLKKNYFVGATPLQKRLELVRRQNSDFSPSPRRKIPQCSQLQEDIDPQKVAFLH KQWTIYSLTPLYKFSYSLKDYSRLLSAFIVAEEKQGVAVEVGEDFNIKVIFSTLLGVKGTQRDHEAFLV QILSKSQSSREHREDKVLWTGWFCVFGESLQETVSEDFCLPLFLANGAESNTSLIRDWFQKTFDCCFS PLAISAFNLSWMAAMWTACKMDRYMATTEFLWSVPCSPQSLDISYAIHPEDAKALWESVHKTPGEVTQE E VDLFMNCLYSHFHRHFKIHLAATRLVRVSTSVASAHTDGKIKILCHKYLIGVLAYMTELAIFQIE TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	39.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001153402
Locus ID:	70454
UniProt ID:	Q14A61


[View online »](#)

RefSeq Size:	2327
Cytogenetics:	1 H2.1
RefSeq ORF:	1035
Synonyms:	2610300B10Rik; AW121806; AW550697; CENP-L