

Product datasheet for TP518239

Zc3h12c (NM_001162921) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse zinc finger CCCH type containing 12C (Zc3h12c), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR218239 representing NM_001162921 Red =Cloning site Green =Tags(s)

MEAGPGAAEGNGPTAVILPLGPSRARGEYGVLCIQEYRKSSKVESSARNFSMGLKDHLGHDLGHLYMEST
 DPQMSAAVPWPMVEKPTMDTVNSGKEGKGVSEENVSSGDSSEGSTSSDHESEQLSSLSVEPCSLTKTHRQ
 L
 CRSPCLEPRLKHS DILQDFKPEESQTPSKEVKKPPDVVREYQTKLEFALKLGYSEEQVQLVLNKLGTDA
 LINDILGELVKLG NKSEADQTVSTINSVMRETSSLESQRSESPMQEVVDDGENLRPVVIDGSNVAMSHG
 NKEVFSCRGIKLAVDWFLERGHKDITVFVPAWRKEQSRPDALITDQEILRKLEKEKILVFTPSRRVQGRR
 VVCYDDRFIVKLAFESDGIIVSNDNYRDLANEKPEWKKFIDERLLMYSFVNDKFMPPDDPLGRHGPSLDN
 FLRKKPIVPEHKKQPCPYGKKCTYGHKCKYHPERGSQPQRSVADELRAMSRNTAAKTTNEGGLVKSNSV
 PCSTKADSTSDV KRGAPKRQSDPSIRTHVYQDIEEKLPTKNKLETRSVPSLV SIPATSTAKPQSTTPLSN
 GLPSGVHFPPQDQRPQGQYPPMMMATKNHGTMPYEQYPKCDSPVDVGYY SMLNAYSNL SISGPRSP
 ERR
 FSLDTDYRVNSVASDCSSEGSMSGSSDSYVGYNDRSYVSSDPQLEESLKCQHMHPSRLNSQPFLQN
 F
 HDPLTRVQSYSHEEPKFHPKRPLPHLAMHLQHPAVGARSSCPGDYPSPPSSAHSKAPHLGRSLVATRIDS
 ISDSRLYDSSPSRQRKPYSRQEGLGSWGRPSYGLEAYGYRQTYSLPDNSTPPCYESITFQSLPEQQEPTW
 RIPPYCGMPHDPPRYQDNREKIFINLCNIFPPDLVRLVMKRNP HMTDAQQLAAAILVEKSQ LGY

TRTRPLEQKLISEEDLANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	101.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001156393</u>
Locus ID:	244871
UniProt ID:	<u>Q5DTV4</u>
RefSeq Size:	6133
Cytogenetics:	9 A5.3
RefSeq ORF:	2709
Synonyms:	A230108E06; C230027N18Rik; mKIAA1726
Summary:	May function as RNase and regulate the levels of target RNA species.[UniProtKB/Swiss-Prot Function]