

Product datasheet for TP518278

Lrrn4cl (NM_001013019) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Purified recombinant protein of Mouse LRRN4 C-terminal like (Lrrn4cl), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug

Species: Mouse

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >MR218278 protein sequence
Red=Cloning site **Green**=Tags(s)

MLGSLSLWLAAMTTSLVSQPQILTLEDYQEGEEDDVTVATPSLAVRCDYDRCRHLQVSCQELQKVG PVA
CLCPGLSREDQQPEPRLGEVQIMAE EGYAVVHWCAPFSPVSHYWLLLWESNGAPQKSAPLNATVRR AEL
KGLKPGVAYVLCVVAANDAGESNVPGA EVEG PENWTGPSFGPCRKFIMPPKPVTLVYAAVGVGTALALLS
CAALVWHFCLRERWGCPRRQGMAQASEAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-MYC/DDK

Predicted MW: 25.9 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001013037](#)

Locus ID: 68852

UniProt ID: [Q3TYX2](#), [Q8C2Z8](#)

RefSeq Size: 2583



[View online »](#)

Cytogenetics: 19 A

RefSeq ORF: 720

Synonyms: 1110067112Rik