

Product datasheet for TP518437

Ccdc159 (NM_025977) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse coiled-coil domain containing 159 (Ccdc159), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA	>MR218437 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MNGCDSLLAGLSRDSDDCVSHLCTPPLLAITSWAAKERSQYSPPPGEPGPSSDKALDCMTHWSRVLEPEI
VNTASEQAGKSGAWKEWDSEPPHEGTPCSSSDVNKDHYHDQDLVKRNHCVAKKSLEPSSAKVKVKNM
IIPDSQKLLRCELSRSQLQAQSKAFEFLNHSVTMLEKESCLQQIKIQLEEVLSPTSROGKEYGRKWS
TEQELYGALAQGLQGLQKTLKEGEELQRARTTRCLQLLAREIRDSKKFLWEELELVREEVTFIYQKLQDQ
EDEISENLLNIQKMQKTQVKCRKVLTKMKQQAYDSWPEAEGVPTEGNGCCKDDLQKELGDIWSAVHSLQS
SIDCLALSMGTRPRASSLRGQKGHQCKSSQCPSWSDSDSDWERPFKSKSGSYPPA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	45.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_080253
Locus ID:	67119
UniProt ID:	Q8C963



[View online »](#)

RefSeq Size: 1531
Cytogenetics: 9 A3
RefSeq ORF: 1212
Synonyms: 2510048L02Rik