

Product datasheet for TP518510

Ccdc157 (NM_001164620) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse coiled-coil domain containing 157 (Ccdc157), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR218510 representing NM_001164620 Red=Cloning site Green=Tags(s) MAHLLGSQACMDSLRLKDLTDLQGTIVDVFSRAGPVRFPSPWKFPDRVACDLDMVALLEHYDHVPGDPEFT Q LSHAVLLELVIDRLLLLLQSCASYLENLSVEQMMPPARAAGPCMSVGLTVRRFWSNLLRLGLLYQQAVPQ KRANQGEISITKPTAKGEPARSPECMTAKFIKPPSPVGLPLICQGLQSIPVRVSLRSPGGTSEKTKSVY SQTVETALVPCDACTSVQGSLWEVGKVVISLCQSQNLPSGLGQFQKLVKDSLGLKPLPAATVGHWAEEQS KDLTRLNKHVGALTQLVGPLRAQLEDAEGQKDGLRKQVSKLEQALQQEQGQRQRQTEEAERTLAKCEHD R HQLLTETCDLKTAVLEGLDKQQQKSIQAMEAKAQQLLEEGERRAAERQVQQLEEQVQLLAGRLDGA S QQIRWASTELDKEKARVDSMVRHQESLQAKQRTLLQQLDCLDQEREELRGSLEAEAAQRSELEEQLQSLQ SDREQEQCQLQAQQELLQSLQKEQDLEQVTTDLQLTISELRQLEELKERERLLVAFPDHLHQPEEAQIQ SSSNVTQDMERQVQANAIQVLQEENKRLQSMILTIREVAQQGGLKMVPQGQLWSPPYKGIQGATPP AQ AQSAFSGLTGRRQSPGSRSTSSTGRTHPGGLRTSPSRQPGGLPSKFSGLGDGSHSASCTQNPICALRLRRK LSPNRQAGSTYQPQERPT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	80.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol


[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001158092</u>
Locus ID:	216516
UniProt ID:	<u>Q5SPX1</u>
RefSeq Size:	5017
Cytogenetics:	11 A1
RefSeq ORF:	2154
Synonyms:	4930562D19Rik; mKIAA1656