

Product datasheet for TP518711

Avl9 (NM_030235) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse AVL9 cell migration associated (Avl9), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR218711 protein sequence Red=Cloning site Green=Tags(s)

MEKSGRES DGAPCGPVLHIVVVG FHHKKG CQVEFSY PPLIPGDG HDSHTLP EEWKYL PFLALPD GAHNY
 Q
 EDTVFFHLPPRNGNGATVYGISCYRQIEAKALKVRQADITRETVQKSVCLSKLPLYGLLQAKLQLITHA
 YFEEDFSQISILKELYEHMNSSLGASLEG SQVYLGLSPRDLVLFHRHKVLILFKLILLEKKVLFYISP
 VNRLVGALMTVLSLFGMIEHGLSDCSQYRPRKSMSEDAGPQESNPSADDFTESTSDVLNTSLETVTRV
 MAVNHGEDAVPKTEKPYFQVEGNNGKQEPSSDGRYLELPPRSPESSESDWETLDPSVLEDASLKEREQ
 MGSDQTHLFQKDSLPSDSPITVQPQANNRQVWLIPGLISGLEEDQYGMPLAIFTKGYLCLPYMALQQHH
 LLSDTVTRGFVAGATNILFRQQKHLSDAIVEVEEALIQIHDPRLKLLNPTTADLRFADYLVRHV TENRD
 DVFLDGTGWEGGDEWIRAQFAVYIHALLAATLQLDNEKMLSDYGTTFVAAWKNTHNYRVWNSNKH PAL
 SE
 INPNHPFQGGQYSVSDMKLRFSHSVQNSERGGKIGSVMVTTSRNVVQTGKAVGQSVGGAFSSAKTAMSS
 WL
 STFTTSTPQSLPEPPNGKP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	72.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


[View online »](#)

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_084511</u>
Locus ID:	78937
UniProt ID:	<u>Q80U56</u>
RefSeq Size:	6691
Cytogenetics:	6 B3
RefSeq ORF:	1947
Synonyms:	5830411G16Rik; D730049P16Rik; mKIAA0241
Summary:	Functions in cell migration.[UniProtKB/Swiss-Prot Function]