

Product datasheet for TP518955

Tsnaxip1 (NM_024445) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse translin-associated factor X (Tsnax) interacting protein 1 (Tsnaxip1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR218955 representing NM_024445 Red =Cloning site Green =Tags(s)

MANLQERKSFSPRISIQASGGTPEAKGIEKRKLSQKRRTLPLQLSLGGHLSWPPTYTSGQTVLHNRKPC
SDDSRNRANSCQQQSMSISKPKYLEQLENYLKELLLDLSTDSAQELRLQPYREIFFEFDKTYKPL
LSSIKNAYEVM LAHQKERIRSLEPLKAKIVTMNEDCSERVLAMRAEERYEISVLKKEKMNLLKLIDKNE
EKISLQSEVAKLRRSLAEEYLRYLTERDARKILIGDLNELRYQQEDMSLAQTPGVWGEDPVKLTALKMT
RQDLTRTQMELNMTKANSGDVVPRRDLEMQEKTNMELQEQLSKADYEEVQKEHELLLQLHMSTLKERD
QFYNELQEIQRTSTPRPDWTKCESIAGGPERWLVAEGKNSDQLVDVLLLEEIGMGLLREKDFFPGLGFG
DAIPPFLRFDGPVKNKKPSKKEVVNLLKDAWKERIAEEQKEPFPDFFNFLERRFGVNDAMAWAYTIFEN
IKLFRSSEIMNQFYAVLMGKNLESVYINQKKTLSHLLKELLSVDTQNEGSITMEQFSTILKTTFFLKKEE
QIQELMEAVGWGPDSSNTDMLNYRSLFNEDEEGQSEPFVQRLWEQYESDKEAYLEELKQELDLDPLEDVT
LLKMRGTLMNIDPTMDKQTL SAYLSQAYQIP AIDVPLEDEEKQGIISIKVETALDRLRMADTKRVGPREP
DPAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	81.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_077765</u>
Locus ID:	72236
UniProt ID:	<u>Q99P25</u>
RefSeq Size:	2210
Cytogenetics:	8 D3
RefSeq ORF:	2112
Synonyms:	1700016K08Rik; TXI1
Summary:	Possible role in spermatogenesis.[UniProtKB/Swiss-Prot Function]