

Product datasheet for TP518986

Zc3hav1l (NM_172467) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse zinc finger CCCH-type, antiviral 1-like (Zc3hav1l), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR218986 representing NM_172467 <div> <div>Red</div>=Cloning site <div>Green</div>=Tags(s) </div> <p> MADLTLC AFLTKVLCAHGGRMFLQDLRGHV LSEAKLR AVLRRAGPERFLLQ EVELRDGPWD AEA EVAAG EGAGAGSGGGATACRVM AVSSARLCARYQRGECQGCDQLHLCRRHMLGKCPHRDCWSTCSLSHDIHTP IN VQVLKSRGLFGLNEGQLRILLQNDPCLFPEVCQMYNKGVDVLYGYCSLKDR CNKFHVCKSFVRGECPFQ PCKRSHQLIHAATLKLEDQELSVSSVNFQIISVYRHKKLHKMLEEKDHSASTEQPQGLGKQGALGAVE ARPFLPARAQSPRPQ </p> <div> <div>TR</div> <div>TRPLEQKLISEEDLAANDILDYKDDDDKV</div> </div>
Tag:	C-MYC/DDK
Predicted MW:	33.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_766055
Locus ID:	209032
UniProt ID:	Q8BFR1


[View online »](#)

RefSeq Size:	6538
Cytogenetics:	6 B1
RefSeq ORF:	888
Synonyms:	B130055L09Rik; E430016P22Rik