

## Product datasheet for TP519084

### Lrrc28 (NM\_175124) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse leucine rich repeat containing 28 (Lrrc28), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR219084 representing NM_175124 Red=Cloning site Green=Tags(s)
	MASEICKTISVARLEKHKNLFLNYRNLHHPLELLKDEGLQHRLERLYMKRNSLTTLPENLAQKLPNLVEL YLHSNNIIVVPEAIGSLVKLQCLDLSDNAIEIVCPEIGGLRALRHLRLANNQLQFLPPEVGDLELQTLTLD ISSNRLALPERLHCLSLQYLTVDNRNLCCVPRHLCQLPSLNELSMAGNHLASLPIDLGRSRELQYVYV DNNIQLKGLPSYLYNKVIGCNGCGIPIQLSEVRLTFFSSGQLTVFLPAEVKTIGTEKDHVLPQLQELTMRS LYRTYHGLWKDLNFLSPISLPRSLLELLHCPLGHCHLCSEPMFTFVYPKIFPLRETPMAGLHQRRTSIGF VAYCCSTQCLRTFNLLC
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	41.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_780333</a>
Locus ID:	67867
UniProt ID:	<a href="#">Q3TX51</a>



[View online »](#)

RefSeq Size: 2978  
Cytogenetics: 7 C  
RefSeq ORF: 1101  
Synonyms: 1300004K21Rik; 2210012C09Rik; 2310058O11Rik