

Product datasheet for **TP519392**

Brox (NM_027861) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse BRO1 domain and CAAX motif containing (Brox), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR219392 representing NM_027861 Red =Cloning site Green =Tags(s)
	MTHWFHRNPLKATAPVSFNYYGMITGPPASKICNDLRSARTRLLELFTDLSCNPETMKNAADLYFSLQGFINSVGDSTQESKLRYIQNFKWTDTLQGHVPSAQQDAVFELISMGFNVALWYTKYASRLAGKENITEDEAKEVHRSLKIAAGIFKHLKESHIPKLLTPAEKGRDLEARLIDAYIIQCQAEAEQVTIARAIELKHAPGLIALAYDTASFYQKADHTLSSLEPAHSKWRKYLHLKMCFYTAYAYCYHGQTLLASDKCGEAIRSLQEAELYAEAEALCKEYGETKGPPTAKPSGHLFFRKLGLSVKNTLDKCQRENGFIYFQKIPTAPQLKANYGLVEPVFPFPPMSAHWTPEALAAFDLTKRPKDDSVKPKPEEDVKPKPEDIRPQKDTGCSVS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	46.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_082137
Locus ID:	71678
UniProt ID:	Q8K2Q7



[View online »](#)

RefSeq Size:	3457
Cytogenetics:	1 H5
RefSeq ORF:	1233
Synonyms:	0610010K06Rik