

Product datasheet for **TP519780**

Fsp1 (NM_027759) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse fibrous sheath-interacting protein 1 (Fsp1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR219780 protein sequence Red =Cloning site Green =Tags(s)
	MPMDIIKGNLDGISKPASSRSRPGSQSSNGSLEVLTPPEGPVKIDMVNKLNSGQEGHTSNNGVEERRNS NDAKWADDSTKPKVKESSEDPDMPQPQATPEHSDDPKLEETNAVLQNAIRKMHRLDKLLAKKQCREK EV KKQGLEMRVKLWEELKSAKNTEDLENDEELGNTKKFLCLTSESAGKAAAEALHCEFEDALFSVFHTQIPP ETYENLMEKDFTCDVEKNEPLITEKQPFSENTEAIEPRSEDSQGFIRQNAEHSQDFIKRNIELAKHSRSP VVMVEGEKKRLDELLQGLDDADSLSSAEGDQCGWLVPGEYTLAATESQQLAEIDIKLQELSVDSPITF SLESQSHKGDMECDANEERNTEPTGKILQDRKEQRDRESRLRAIDGKLKELSEQVEECPMITPSKRNE RITWRWLLAKILEPE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	49 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_082035</u>


[View online »](#)

Locus ID: 71313

UniProt ID: [Q9D3V5](#)

RefSeq Size: 1673

Cytogenetics: 2 E5

RefSeq ORF: 1305

Synonyms: 1700012M13Rik; 4933432K11Rik