

Product datasheet for TP519952

Sash3 (NM_028773) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse SAM and SH3 domain containing 3 (Sash3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR219952 representing NM_028773 Red =Cloning site Green =Tags(s) MLRRKPSNASDKEPTQKKKLSLQRSSSFKDFAKSKPSSPVVSEKEFNLDDNIPEDDSGVLTPEDSGKSGK KLGKKWRAVISRTMNRKMGKMMVKALSEEMGDTLEEGSASPTSPDCSLDSPGPEKMALAFTEQEEREPP S LSRQTSTGSELCSPPGSGSFLEESPAPQYTGPFGRARVHTDFTSPYDHDLSLKLQKGDVIQIVEKPPV GTWLGLLNGKLGSFKFIYVDVLP EEAVGPVRPSRRQSKGKRPKPTLHELLERIGLEEHTSTLLNGYQT LEDFKELRETHLNLNIMDPQHRAKLLTAAELLDYDTGSEEAEEGAESSQEPVAHTVSEPKVDIPRDSG CFEGSESGRDEAELAGTEQLQGLSLSGAP TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	41.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_083049
Locus ID:	74131


[View online »](#)

UniProt ID: [Q8K352](#)

RefSeq Size: 2692

Cytogenetics: X A4

RefSeq ORF: 1140

Synonyms: 1200013B08Rik; AW413946; SLY1

Summary: May function as a signaling adapter protein in lymphocytes.[UniProtKB/Swiss-Prot Function]