

Product datasheet for TP520199

Lptm4a (NM_008640) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse lysosomal-associated protein transmembrane 4A (Lptm4a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR220199 representing NM_008640 Red =Cloning site Green =Tags(s)
	MWQKAEGEIPSRPARAAPGAQSAGAQLARARAEGGPEDASFGCGHRRRVWVVEPLGRKLLSGAWARTR PARFRMVSMFKRSRDRFYSTRCCGCCHVRTGTIILGTWYMVVNLMAILLTVETHPNMSPAVNIQYE VIGNYSSERMADNACVLFAVSVLMFISSMLVYGAIYQVGWLIPFFCYRLDFVLSCLVAISSLTYP RIKEYLDQLPDFPYKDDLALDSSCLLFIVLVFFVVFIIFKAYLINCWVNCYKYINNRNVPEIAVYPAFE TPPQYVLPTYEMAVKIPEKEPPPPYLPA
	SGP TRTRPLE QKLISEEDLAANDILDYKDDDDK V
Tag:	C-MYC/DDK
Predicted MW:	35.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_032666
Locus ID:	17775
UniProt ID:	Q8BG66 , E9QLR3



[View online »](#)

RefSeq Size: 2130
Cytogenetics: 12 A1.1
RefSeq ORF: 924
Synonyms: AA286466; LAPTM4; mKIAA0108; MTP; Mtrp