

Product datasheet for TP520390

Fstl4 (NM_177059) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse follistatin-like 4 (Fstl4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR220390 protein sequence Red =Cloning site Green =Tags(s) MKPGGFWPHLALLGVSLPAVLGWMDQGASRSPNMVPGESQAETRGRFEVTRREGLPSPGLSASCGKKL CS HGSRCLLNRTTGQPSCQCLEVCPRYPMPVCGSDGRLYGNHCELRRACLKGKRVSVHSKDCFLKGDMC T MAGYARLKNVLLALQSRRQPLPQGTPRQSLASQKRLLVESLFKDLADGNGHLGSLELAQYVLKEQDMD G SLNRCSPDLLRFDDYNSDGSLLTGEFYTAQVIQLSLAPEDKVSVTTVTVGLSTVLTCAIRGDLRPPII WKRNGLTLSFLGLEIDINDFGEDGSLYITKVTTVHMGNYTCHALGHEQLVQTHVLQVNVPPVIRVYPETQA QEPGVAASLRCHAEGIPLPRIIWLKNGMDVSTQMSKQLSLLANGSELHIGSVRYEDTGAYTCIAKNEGV DEDISSLFIEDSARKTLANILWREEGLSVGNMFYVFAEDGIVVIHPVDCEVQRRLKPTKIFMSYEEICP RVEGDATQPCQWASAVNVRNRYIYAQPALNRVLVDVQAQKVLQSIGVDPLPVKLSYDKSHDQVWVLS W GDMHKSQPSLQVITEASTGQGQHLIRTPFAGVDNFFIPPTNLIINHIFGFIFNKTDPAVHKVDLETMA LKTISLRHYGCMPQAMAHTLGGYFFVQCQQDTPSTGTPQLLIDSVTDSVLGPNSDITGTPHVSPDGRFI VSVSNKGPWLHVQEVTVRGEIQTLYDLKINPGISDLAFQHSFTEGSQYNAYATLDKEPDLLFLELSTGKM GRLKNLKEPPRGPAPTWGGPRRVLRDSEGLFGQYLLTPAQESLFLVNGRQNALRCEVSGIKGAATVWWGE V TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	92.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining


[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_796033</u>
Locus ID:	320027
UniProt ID:	<u>Q5STE3</u>
RefSeq Size:	3117
Cytogenetics:	11 B1.3
RefSeq ORF:	2523
Synonyms:	B230374F23Rik; mKIAA1061; SPIG1