

Product datasheet for TP520446

Sowahd (NM_173779) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse sosondowah ankyrin repeat domain family member D (Sowahd), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR220446 representing NM_173779 Red =Cloning site Green =Tags(s)

MAQALEDGNPLPKASNRPAESEAPSDPQIKDNHCLGRYKQAVRDSANLSQERIFHSALTVSGASWARRR
GELRELLGLQGAAPAGWLSEEHLEPAVPGSRRSSGQSSRVCLEPREHAWILAAACRFEVLLLEMLEAEP
SLLMREDPITGYSVLHWLAKHGRHEELILLHDFARRRGLPFDVSAPGSGGLTPLHLAALQGHDMMVIKVLV
GALGADPSRRDHSGNRPCHYLRPDASLNLRELSGAEWEIERDRKRENANNSSRTTTTTTTTSRWLKR
PSASCIKSTGVHYKEASQPVKEKKASSTQEGQGRCLRRYLPFTQNR

SGPTRTRPLE**QKLISEEDLAANDILDYKDDDDK**V

Tag:	C-MYC/DDK
Predicted MW:	36.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_776140
Locus ID:	245381
UniProt ID:	Q8BY98



[View online »](#)

RefSeq Size: 1574
Cytogenetics: X A3.3
RefSeq ORF: 981
Synonyms: A630014H24Rik; Ankrd58