

Product datasheet for TP520693

Klhdc9 (NM_001033039) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse kelch domain containing 9 (Klhdc9), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA	>MR220693 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MAGVQPLGRARGSTWTWRPVARDVLLARAFHSCTELEGRFYLVGGLLEGGARVPSNDTVIFDPAVGQAVR
LVARGSPLRSHHDAALVGGRWLCVGGWDGSRRLSTVAALDTERGVWEAWAANPGNCPAGLSSHTCTRL
SDGELRVSGREGGHTQRRYGSYITLKLDRHRTCYKEEGCHTTSRSGHCAALLPTAGPHPGHQLLLFG
GCNSVGPVAGQWSPGKIKEEQPVAPHLREQLARLVSSGQGLQGPQSLRHHSCSVGPFVAVLFGGETLT
RARDTICNDLYIDTRKSPPLWFHFPSTDRGLKRVGHRTCLWNDQLYLVGGFGEDGRTASPQVCILEFFI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	38.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001028211
Locus ID:	68874
UniProt ID:	Q3USL1



[View online »](#)

RefSeq Size: 1555
Cytogenetics: 1 79.36 cM
RefSeq ORF: 1053
Synonyms: 1190002J23Rik; AA087357; ESTM31