

Product datasheet for TP520728

Spata2l (NM_030176) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse spermatogenesis associated 2-like (Spata2l), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR220728 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGSSSLSEDYRQCLERELRRGRAGVCGDPSLR AVLWQILVEDFDLHGALQDDALALFTDGLWGRADLAPA LQDLARAFELLELA AVHLYLLPWRKEFTTIKTFSGGYVHVLKGV LSEELLTRSFQKMGYVRRDNHRLMVT TPPPACQLVQVALGCFALRLECEILSEVLTQLGTSVLP AEELLRARRASGDVASCVAWLQQR LAQDEEPP PLPPRGTPATYGAPVDLYQDLQEDESSEASLYGEPSPGLD SPPVELAYRPPLWEQS AKLWGS GGQPWEPP ADDMHRASSPPYGALEEELEPEPSAFSFLSLRRELSRSGDL APPESSPGQASPRHRQAEAAAASSAYGP AVEPLSYQAHSCLSPGNLPTLCCDTCRQLHATHCTALSACRPTHSLRILLGDNQRRLWLQRAQVDNLLYD SPGAHP</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-MYC/DDK
Predicted MW:	46.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_084452</u>
Locus ID:	78779



[View online >](#)

UniProt ID: [Q8BNN1](#), [B2RTE8](#)

RefSeq Size: 2380

Cytogenetics: 8 E1

RefSeq ORF: 1281

Synonyms: 2610039E05Rik; AI648233; mKIAA4138