

Product datasheet for TP521105

Rnf157 (NM_027258) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse ring finger protein 157 (Rnf157), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221105 representing NM_027258 Red =Cloning site Green =Tags(s) MGALTSRQHAGVEEVDIPSNSVYRYPKSGSYFASHFIMGGEKFDSTHPEGYLFGENSNDLNLGNRPVAF PYAAPPPQEPVKTLRSLINIRKDTLRLVKCAEEVKSHGEEAGKAKVHYNVEFTFDTDARVAITIIYYQATE EFQNGIASYIPKDNSLQSETVHYKRGVCQQFCLPSHTVDPSEWAEELGFDLDREVYPLVVHAVVDEGDE YFGHCHVLLGTFEKHPDGTFCVKPLKQKQVVDGVSLLQEIYGIENKYNTQDSKVAEDDVSDNSAECVVC LSDVRDTLILPCRHLCLCNTCADTLRYQANNCPICRLPFRALLQIRAMRKKLGPLSPSSFNPIISSQTS SEEHSSSENIPPGYEVVSLLEALNGPLTSSPAVPPLHVLGDGHLGMLPSYGSDDGYLPPVRTLSPLDRLS DCNNQGLKLLKSLSKSISQNSSVLHEEEDERSCESDTQLSQRLSVQHPPEEGPDVTPESENLTSSSGAV DQSSCTGTPLSSTISSPEDPASSSLAQSVMSMASSQISTDTVSSMSGSYIAPGTEEEGEALPSPRAASRA PSEGEETPAESPDSNFAGLPAGEQDAEGNDIIIEEDRSPVREDGQRTCAFLGMECDNNNDFDVASVKALD NKLCEVCLPGTWQHEDNTVNCRHTQRRRLSSSSLEDPEENRPCVWDPMVA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	74.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: [NP_081534](#)
Locus ID: 217340
UniProt ID: [A2AAN9](#)
RefSeq Size: 4655
Cytogenetics: 11 E2
RefSeq ORF: 2043
Synonyms: 2610036E23Rik; A130073L17Rik; mKIAA1917